

## RELY™ H. pylori RAPID TEST

After a brief and transient clinical illness, H. pylori infections are asymptomatic, yet clinical trials have demonstrated that H. pylori is associated with 90% of patients with duodenal ulcer and approximately 80% of those with gastric ulcer.<sup>1</sup>

Help reduce the risk of a more serious illness with...

### Fast, accurate and reliable detection of H. pylori

<sup>1</sup>CDC, Helicobacter pylori Fact Sheet for Health Care Providers.



### Easy and reliable detection of H.pylori means enhanced patient outcomes...

- Minimizes the use of inappropriate treatment protocols
- Early detection and effective management of H. pylori reduce patient visits and cost of clinical follow-up
- Helps lower the risk of developing a more serious illness

### RELY™ H. pylori Rapid Test provides simple, fast and cost-effective detection of H. pylori.

- **Simple and Easy-to-use** – CLIA waived\*, requires minimal training
- **Rapid Results** – Two step test provides results in 10 minutes
- **Flexible** – Multiple sample types (Whole Blood, Serum, or Plasma)
- **Reliable** – Procedural and external controls provided
- **Cost Effective** – Convenient room temperature/refrigerated storage and extended dating help manage inventory during seasonal shifts in testing

\*Whole blood

\*\*Whole blood, Fingerstick

<b>Format</b> .....	20 Cassette Test Kit
<b>CLIA Status</b> .....	Waived: Whole Blood Non-Waived: Serum and Plasma
<b>Sample</b> .....	Whole Blood (Fingerstick or Venipuncture), Serum or Plasma
<b>Sensitivity**</b> .....	.87%
<b>Specificity**</b> .....	.86%
<b>Procedural Steps</b> .....	2
<b>Time to Results</b> .....	10 Minutes
<b>Shelf Life</b> .....	18 Months from Date of Manufacture
<b>Storage</b> .....	2-30° C (36-86° F)



Product Description	Size	Reference No.	CPT Code
RELY™ H. pylori Rapid Test	20 Test	6300-020	86318QW (Whole Blood) 86318 (Serum/Plasma)



- General & Special Chemistry
- Pregnancy
- Point-of-Care
- Rapid Tests
- Therapeutic Drug Monitoring
- Serology
- Occult Blood
- Urinalysis
- Clinical Instrumentation
- Pipettors

1261 North Main Street  
Boerne, Texas 78006 USA  
Toll Free: 866.STANBIO  
(782.6246)  
Tel: 830.249.0772  
Fax: 830.249.0851

stanbio@stanbio.com  
www.stanbio.com



DISTRIBUTED BY