

QuStick™ hCG

Pregnancy Test Strips

Sensitivity: 20 mIU/mL in Urine



The Stanbio Laboratory QuStick™ Pregnancy test is a qualitative immunoassay for the detection of human chorionic gonadotropin (hCG) in urine.

- CLIA Waived
- Long shelf-life
- Built-in (+) & (-) procedural controls
- Accuracy: > 99%
- Urine specimens
- Dip & read strip format
- Positive results in as little as 30 seconds
- Specificity: > 99%

Ref. No.	Size
1210-425	100 Strips (4 re-sealable canisters, 25 each)

Stanbio Laboratory • 1261 North Main Street • Boerne, Texas 78006
Toll Free: 866-782-6246 • Tel: 830-249-0772 • Fax: 830-249-0851
Web: www.stanbio.com

Distributed by



RSS0023.00