

Q U S T I C K S T R E P A



Ultra-sensitive.

Stanbio's QuStick Strep A detects many tough cases others miss since fewer colonies are required to yield a positive result.

- 91% @ rare (< 5 colonies)
- 100% @ 1+ (< 10 colonies)

Convenient. Economical.

QuStick eliminates waste with room temperature storage and a two-year shelf life from date of manufacture. You get results in five minutes and, QuStick is not only priced less, it outperforms the leading brand.

QuStick Strep A

Each kit includes

- Individually pouched test strips
- Reagent A and Reagent B
- Positive and negative external controls
- Dacron throat swabs
- Extraction tubes
- Instructions for use
- Handy pictorial quick guide showing step-by-step procedure.



Reference Nos. 6000-025 (25 Tests)
6000-050 (50 Tests)



CODE*: 87880QW

*Exemplar code provided. Refer to regional CMS fee schedule for relevant codes and reimbursement rates.



STANBIO
LABORATORY

The World's Source For Clinical Diagnostics

- General & Special Chemistry
- Pregnancy
- Point-of-Care
- Rapid Tests
- Therapeutic Drug Monitoring
- Serology
- Occult Blood
- Urinalysis
- Clinical Instrumentation

Stanbio Laboratory
1261 North Main Street
Boerne, Texas 78006 USA

866.782.6246 (U.S. Toll Free)
Tel: 830.249.0772
Fax: 830.249.0851

email: stanbio@stanbio.com
www.stanbio.com



DISTRIBUTED BY

RSS0008.03
Eff. Date: 4/19/11
DAR #: 11.0272
©2011 Stanbio Laboratory
Printed in USA.